

MADA EMERGENCY OXYGEN KITS

No. 1300ME




No. 1301BME



No. 1301ME



No. 1528E



No. 1530E



No. 1528BE



No 1629E



No. 1631E



Mada Medical Products, Inc.

DOMESTIC / INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 Fax: 201-460-3509

To Order Call Toll Free: 1-800-526-6370

To order unit filled, remove suffix "E"

Visit www.madamedical.com for a complete listing of our products.

MADA EMERGENCY OXYGEN UNITS

PRODUCT NO.	CYLINDER	REGULATOR	CASE/CART	FEATURES
1300ME	(M7) 1303ME	1308A 6 LPM Avg.	1319A	Inhalator only
1301ME	(M7) 1303ME	1308A 6 LPM Avg.	1318A	Inhalator only
1301BME	(M7) 1303ME	1308A 6 LPM Avg.	1410	With BMR
1501E	(MD) 1503E	1442 6 LPM Avg.	1520	Inhalator only
1528E	(MD) 1502E	1441 6 LPM Avg.	1500A	Inhalator only
1528BE	(MD) 1502E	1441 6 LPM Avg.	1500A	Inhalator with BMR
1529E	(MD) 1502E	R1835-25GB (1/4-25 LPM)	1520	Inhalator only
1530E	(MD) 1502E	R1835GBHF-25 (1/4-25 LPM)	1500A	Demand Valve with Inhalator
1629E	(ME) 1603E	1442 6 LPM Avg.	370C	Inhalator with BMR
1629DE	(ME) 1603E	1442-10 10 LPM Avg.	370C	Adult/Child BMR
1630ADE	(ME) 1602E	1335 (2-8 LPM)	370C	With BMR
1631E	(ME) 1602E	1635-1 (50 PSI Reducer)	370C	Demand Valve with no Inhalator
1634E	(ME) 1602E	R1835GBHF-25 (1/4-25 LPM)	370C	Demand Valve with Inhalator

ACCESSORIES

1432	Reusable Adult bag mask resuscitator (Each)
4000	Disposable Adult bag mask resuscitator (Box 6)

To order unit filled, remove suffix "E"

Other valve combinations available upon request.

Filled Kits available upon request.

Note: BMR is abbreviation for Bag Mask Resuscitator used for manual resuscitation.

Demand Valve is used for Oxygen powered Resuscitation.

Inhalator is for administering low flow oxygen only.

Distributed By:



Mada Medical Products, Inc.
DOMESTIC / INTERNATIONAL
625 Washington Avenue, Carlstadt, NJ 07072
Tel: 201-460-0454 Fax: 201-460-3509
To Order Call Toll Free: 1-800-526-6370