

# No. 1701 Series

## MEDICAL FLOWMETERS



PRECISE ADMINISTRATION IN FRACTIONAL INCREMENTS

ADAPTABLE FOR ALL 1/8 NPT INLET ADAPTERS

AVAILABLE FOR OXYGEN OR MEDICAL AIR


5-YEAR WARRANTY

- ✓ DURABLE CYLINDRICAL CONSTRUCTION
- ✓ BACK PRESSURE COMPENSATED
- ✓ POLYCARBONATE HOUSING
- ✓ FLUSH FLOW
- ✓ INLET FILTER

PART NO.	GAS	FLOW RATE (LPM)	ACCURACY (OF FULL SCALE)
1701	O <sub>2</sub>	( 1 – 15 )	± 10%
1701-1	O <sub>2</sub>	( 1/10 – 1 )	± 15%
1701-8	O <sub>2</sub>	( ½ – 8 )	± 10%
1701-A	AIR	( 1 – 15 )	± 10%

# Mada Inlet Adapters

Mada provides a complete line of inlet adapters for the most popular oxygen and medical air flowmeters; factory assembled to your flowmeter upon request.



## ORDERING INFORMATION:

### Product No.

### Description

1701-1706	Oxygen Flowmeter with Ohio Male Inlet
1701-1708	Oxygen Flowmeter with Chemtron Inlet
1701-1708Y	Oxygen Flowmeter with Chemtron Inlet, Y-Branch
1701A-1708A	Air Flowmeter with Chemtron Air Inlet