

Voyage-II Compressor / Nebulizer

from
Mada! No. 277

Delivers effective and efficient aerosol medication therapy with clinical accuracy.

Built for quality & durability for years of reliable treatment.



STANDARD FEATURES

- Powerful piston pump.
- Lightweight: 4 LBS. 2 OZ.
- Storage for accessories.
- Built in handle.
- Maximum pressure: 35 PSI.
- Operating pressure: 18 PSI.
- Liter flow: 6-7 LPM at 11 PSIG.
- Comes complete with nebulizer, tee adapter, mouthpiece and 7' tubing.
- 3 year warranty.
- Electrical requirements: 115 volts, 60Hz

ACCESSORIES

- | | |
|-------|---|
| 1312 | Up - Mist™ Nebulizer: Aerosol "T" with 7' oxygen tubing, single patient use (Box of 50) |
| 50310 | Replacement Medication Nebulizer with adult mask and 7' tubing (Box of 50) |
| 50311 | Replacement Medication Nebulizer with pediatric mask and 7' tubing (Box of 50) |
| 277AF | Replacement Air Filters (Box of 5) |



625 WASHINGTON AVE, CARLSTADT, NJ 07072

1-800-526-6370 • 201-460-0454 • Fax: 201-460-3509

www.madamedical.com