

PENGUIN

Compressor / Nebulizer

from
Mada! No. 278

*The Penguin is a child-friendly
Compressor/Nebulizer that
provides quieter treatment
and durability*



No. 278SB

Optional:
Padded Shoulder Bag

STANDARD FEATURES

- Powerful piston pump.
- Lightweight: 4 LBS. 2 OZ.
- Child friendly.
- Built in handle.
- Maximum pressure: 35 PSI.
- Operating pressure: 18 PSI.
- Liter flow: 6-7 LPM at 11 PSIG.
- Comes complete with nebulizer, mouthpiece, 7' tubing, and mask.
- 3 year warranty.
- Electrical requirements: 115 volts, 60Hz

ACCESSORIES

1312	Up - Mist™ Nebulizer: Aerosol "T" with 7' oxygen tubing, single patient use (Box of 50)
50311	Replacement Medication Nebulizer with pediatric mask and 7' tubing (Box of 50)
277AF	Replacement Air Filters (Box of 5)
278SB	Padded shoulder bag with strap



MADA

MEDICAL PRODUCTS, INC.

625 WASHINGTON AVE, CARLSTADT, NJ 07072

1-800-526-6370 • 201-460-0454 • Fax: 201-460-3509

www.madamedical.com