

# MADA IO<sub>2</sub>

VALVE INTEGRATED PRESSURE REGULATOR



Easy-to-Carry Handle

Easy-to-Read Cylinder  
Contents Gauge

Flow Adjustment Dial (0-25 LPM)

Anodized Barb for Cannula Delivery

Plated Brass Bod

50 PSI Constant Flow DISS Outlet

CGA 540 Fill Port

Portable

- MRI Conditional 3 Tesla
- FDA 510(K) Cleared

## ORDER INFORMATION

Cylinders with the MADAIO2 installed

| ITEM #    | SIZE       |
|-----------|------------|
| 1506IVO2E | MD         |
| 1606IVO2E | ME         |
| MADAIO2   | Valve Only |

## MADAIO2 TECHNICAL INFORMATION

| Items # | Operating PSIG | Pressure | Outlet Flows | Fill Connection | Auxillary Port    |
|---------|----------------|----------|--------------|-----------------|-------------------|
| MADAIO2 | 0-2300 PSIG    | 50 PSIG  | 0-25 LPM     | CGA-540         | DISS 1240 Hi Flow |



1-800-526-6370  
625 Washington Ave.  
Carlstadt, NJ 07072 USA  
[www.madamtedical.com](http://www.madamtedical.com)

